


Drawings


90020 to 90033 TC-K Type Temperature Probes

90020


Unit : mm(approx.inch)



90021




90022




No.	Name	Description	No.	Name	Description
1	Sheath		5		
2	Grip		6		
3	Cable	Heat resistant PVC			
4	Connector				

90023




90024




No.	Name	Description	No.	Name	Description
1	Sheath		5	Protective cap	PVC
2	Grip		6		
3	Cable				
4	Connector				


90030




90032



90031



90033



No.	Name	Description	No.	Name	Description
1	Sheath		5	Protective cap	
2	Grip		6	Pipe	
3	Cable	Silicone rubber	7	Guard	
4	Connector		8	Temperature sensor	