


Drawings

CL150/CL155 Clamp-on Tester

Unit: mm
(approx. inch)



Weight: approx. 470 g
(including batteries)

Note:

This figure shows the CL150, but there is no difference in exterior from the CL155.

If not specified, the tolerance is $\pm 3\%$. However, in cases of less than 10 mm, the tolerance is ± 0.3 mm.

YOKOGAWA ◆

All Rights Reserved, Copyright © 2016 Yokogawa Meters & Instruments Corporation
Printed in Japan

SD CL150-E
4th Edition : Apr. 2016 (YMI)