


Drawings

Model 2792A (2792A01 to 2792A08) Standard Resistor

Unit: mm
(approx. inch)



Weight: approx. 1.2 kg

Model	Nominal Value
2792A01	0.001 Ω
2792A02	0.01 Ω
2792A03	0.1 Ω
2792A04	1 Ω
2792A05	10 Ω
2792A06	100 Ω
2792A07	1 k Ω
2792A08	10 k Ω

If not specified, the tolerance is $\pm 3\%$. However, in cases of less than 10 mm, the tolerance is ± 0.3 mm.