



**90020B, 90021B, 90022B, 90023B, 90024B,
90030B, 90031B, 90032B, and 90033B
TC-K Probe**

90023B




90024B




No.	Name	Description	No.	Name	Description
1	Sheath		5	Protection cap	PVC
2	Grip				
3	Cable				
4	Connector				


90030B




90031B



90032B



90033B



No.	Name	Description	No.	Name	Description
1	- - -		5	Protection cap	
2	Grip		6	Pipe	
3	Cable	Silicone covered cable	7	Guard	
4	Connector		8	Temperature sensor	