


# Drawings

## Model 53005, 53006 EMISSION THERMOMETER

Unit: mm



Weight: approx. 250 g  
(including batteries)

If not specified, the tolerance is  $\pm 3\%$ . However, in cases of less than 10 mm, the tolerance is  $\pm 0.3$  mm.

**YOKOGAWA** ◆

All Rights Reserved. Copyright © 2005 Yokogawa Test & Measurement Corporation  
Printed in Japan

**SD 53005-E**  
3rd Edition: Aug. 2024 (YMI)