

Instrument Discounts & Special Offers

Yokogawa is happy to provide discounts on a limited selection of instruments. Additional options may be available based on remaining supply. Please note these offers are valid only in the United States. To learn more or enquire about a product, please reach out to us via our [Contact page](#).

Electrical Test Equipment						
Model Code	Qty.	List Unit Price	Special Offer Price	Description	Condition	Discontinuation Status
701250	1	\$ 1,265.00	\$ 750.00	High-Speed 10 MS/s; 12-Bit; Isolation Module for ScopeCorder Series	New	Discontinued; Sept 29th, 2017;
701920	2	\$ 1,700.00	\$ 1,100.00	Differential Probe 12V / 500 MHz	New	Discontinued; Nov 27th, 2020;
702912	1	\$ 440.00	\$ 300.00	Logic Probe	New	Not Discontinued
720220	3	\$ 2,520.00	\$ 1,500.00	16-CH Voltage Input Module	New	Not Discontinued
701267	9	\$ 1,265.00	\$ 750.00	Hi-voltage 100kSps16bit Isolation Module	New	Discontinued; Sept 29th, 2017;
701265	4	\$ 1,640.00	\$ 330.00	500S/s; 16-bit High-Precision Volt/Temp Module (2-ch)	Used (Certified Pre-owned)	Not Discontinued
760901	7	\$ 4,540.00	\$ 3,405.00	30A High Accuracy Element	Gently Used (Certified Pre-owned)	Not Discontinued
760303-03-MV-D/G6/DT/DA/C7/C5	39	NA	\$ 5,000.00	Older WT3000. These are 3 element units with 30A Modules This offer is for Universities only.	Gently Used (Certified Pre-owned)	Discontinued; March 31st, 2021;
WT332E-C1-D/C7/G5/US	2	\$ 5,850.00	\$ 5,440.00	WT332E. This ia a 2 element unit.	New	Not Discontinued
760203-D/C1/C7/EX3/G5/DT/FQ/V1 (WT500)	2	\$ 15,360.00	\$ 7,680.00	WT500 Mid-Range 3-ch Power Analyzer	Gently Used (Certified Pre-owned)	Not Discontinued
DLM3054-D-HE-P4/C1/C8/F01/F02/F03 /F04/F05/G02/G03	1	\$ 24,650.00	\$ 12,325.00	DLM3000 Mixed Signal Oscilloscope; 8-bit; 4 Channels; 500MHz	Gently Used (Certified Pre-owned)	Not Discontinued
DLM3054-D-HE-P4/C1/F01/F02/F03/F04 /F05/G02/G03	1	\$ 23,640.00	\$ 11,820.00	DLM3000 Mixed Signal Oscilloscope; 8-bit; 4 Channels; 500MHz	Gently Used (Certified Pre-owned)	Not Discontinued
WT500 - 760203-D/C1/C7/EX3/G5/DT/FQ/V1	2	\$ 15,360.00	\$ 7,680.00	WT500 Mid-Range Power Analyzer	Gently Used (Certified Pre-owned)	Not Discontinued
WT1803E-5A0-50A3-HE-D/EX3	2	\$ 16,350.00	\$ 13,897.50	WT1800E Power Analyzer with 3-50A elements (options available to add at 15% discount)	New	Discontinued; Sept 30th, 2025;
WT1803E-5A3-50A0-HE-D	1	\$ 15,580.00	\$ 13,243.00	WT1800E Power Analyzer with 3-5A elements (options available to add at 15% discount)	New	Discontinued; Sept 30th, 2025;
WT1804E-5A4-50A0-HE-D	2	\$ 17,990.00	\$ 15,291.50	WT1800E Power Analyzer with 4-5A elements (options available to add at 15% discount)	New	Discontinued; Sept 30th, 2025;
WT1804E-5A4-50A0-HE-D/EX4	1	\$ 19,010.00	\$ 16,158.50	WT1800E Power Analyzer with 4-5A elements (options available to add at 15% discount)	New	Discontinued; Sept 30th, 2025;
WT1806E-5A6-50A0-HE-D/EX6	1	\$ 24,340.00	\$ 20,689.00	WT1800E Power Analyzer with 6-5A elements (options available to add at 15% discount)	New	Discontinued; Sept 30th, 2025;

Optical Instruments						
Model Code	Qty.	List Unit Price	Special Offer Price	Description	Condition	Condition
AQ1100A-HE-D-HPM-USC/VLS	1	\$ 4,300.00	\$ 2,750.00	AQ1100 OLTS MULTI FIELD TESTER; SM1310/1550nm	New	Not Discontinued
AQ2781V-SCC	1	\$ 1,200.00	\$ 600.00	AQ2781V High Power OPM & VLS Module	New	Not Discontinued
AQ7280-HE	1	\$ 2,520.00	\$ 1,550.00	OTDR Main Frame	New	Not Discontinued
AQ2200-331	1	\$ 6,275.00	\$ 1,255.00	Optical Attenuator Module for AQ2200 Series	Used (Certified Pre-owned)	Discontinued; Mar 31st, 2017;
AQ2200132-CC-T6-PA-FCA	2	\$ 11,700.80	\$ 9,945.68	Grid TLS Module (2-Channel)	New	Not Discontinued

Portable Test Instruments						
Model Code	Qty.	List Unit Price	Special Offer Price	Description	Condition	Condition
CA150/TE	4	\$ 2,244.00	\$ 1,000.00	HANDY CAL	New	Discontinued; Sept 30th, 2020;
CL135	1	\$ 252.00	\$ 126.00	Clamp-on Tester (True RMS, ACA, ACV, 600A)	New	Discontinued; Oct 31th, 2019;