

Instrument Discounts & Special Offers

Yokogawa is happy to provide discounts on a limited selection of instruments. Additional options may be available based on remaining supply. Please note these offers are valid only in the United States. To learn more or enquire about a product, please reach out to us via our [Contact page](#).

Electrical Test Equipment						
Model Code	Qty.	List Unit Price	Special Offer Price	Description	Condition	Discontinuation Status
701250	1	\$ 1,265.00	\$ 750.00	High-Speed 10 MS/s; 12-Bit; Isolation Module for ScopeCorder Series	New	Discontinued; Sept 29th, 2017;
701250	11	\$ 1,265.00	\$ 260.00	High-Speed 10 MS/s; 12-Bit; Isolation Module for ScopeCorder Series	Used (Certified Pre-owned)	Discontinued; Sept 29th, 2017;
701920	2	\$ 1,700.00	\$ 1,100.00	Differential Probe 12V / 500 MHz	New	Discontinued; Nov 27th, 2020;
702912	1	\$ 440.00	\$ 300.00	Logic Probe	New	Not Discontinued
720220	4	\$ 2,520.00	\$ 1,500.00	16-CH Voltage Input Module	New	Not Discontinued
701267	13	\$ 1,265.00	\$ 750.00	Hi-voltage 100kSps16bit Isolation Module	New	Discontinued; Sept 29th, 2017;
701265	15	\$ 1,640.00	\$ 330.00	500S/s; 16-bit High-Precision Volt/Temp Module (2-ch)	Used (Certified Pre-owned)	Not Discontinued
701989	5	\$ 1,260.00	\$ 250.00	Logic Probe; 250 MHz / 8-bit	Used (Certified Pre-owned)	Discontinued; Jun 30th, 2012;
701981	12	\$ 870.00	\$ 175.00	Logic Probe; 250 MHz / 8-bit	Used (Certified Pre-owned)	Discontinued; Jun 30th, 2012;
700939	5	\$ 850.00	\$ 200.00	Active FET Probe 10V / 900 MHz	Used (Certified Pre-owned)	Not Discontinued
760901	7	\$ 4,540.00	\$ 3,405.00	30A High Accuracy Element	Gently Used (Certified Pre-owned)	Not Discontinued
760303-03-MV-D/G6/DT/DA/C7/C5	40	NA	\$ 5,000.00	Older WT3000. These are 3 element units with 30A Modules This offer is for Universities only.	Gently Used (Certified Pre-owned)	Discontinued; March 31st, 2021;
WT332E-C1-D/C7/G5/US	1	\$ 5,850.00	\$ 4,970.00	WT332E. This is a 2 element unit. Gently used.	Gently Used (Certified Pre-owned)	Not Discontinued
WT332E-C1-D/C7/G5/US	2	\$ 5,850.00	\$ 5,440.00	WT332E. This is a 2 element unit.	New	Not Discontinued
PX8000-D-HE/B5/C20/G5/M2/P4	2	\$ 43,480.00	\$ 21,740.00	PX8000 Precision Power Scope with 4 Current Modules and 4 Voltage Modules	Gently Used (Certified Pre-owned)	Not Discontinued
760203-D/C1/C7/EX3/G5/DT/FQ/V1 (WT500)	2	\$ 15,360.00	\$ 7,680.00	WT500 Mid-Range 3-ch Power Analyzer	Gently Used (Certified Pre-owned)	Not Discontinued
DLM3054-D-HE-P4/C1/C8/F01/F02/F03 /F04/F05/G02/G03	1	\$ 24,650.00	\$ 12,325.00	DLM3000 Mixed Signal Oscilloscope; 8-bit; 4 Channels; 500MHz	Gently Used (Certified Pre-owned)	Not Discontinued
DLM3054-D-HE/P4/C1/F01/F02/F03/F04 /F05/G02/G03	2	\$ 23,640.00	\$ 11,820.00	DLM3000 Mixed Signal Oscilloscope; 8-bit; 4 Channels; 500MHz	Gently Used (Certified Pre-owned)	Not Discontinued
WT500 - 760203-D/C1/C7/EX3/G5/DT/FQ/V1	2	\$ 15,360.00	\$ 7,680.00	WT500 Mid-Range Power Analyzer	Gently Used (Certified Pre-owned)	Not Discontinued
Optical Instruments						
Model Code	Qty.	List Unit Price	Special Offer Price	Description	Condition	Condition
AQ1100A-HE-D-HPM-USC/VLS	1	\$ 4,300.00	\$2,750.00	AQ1100 OLTS MULTI FIELD TESTER; SM1310/1550nm	New	Not Discontinued
AQ2781V-SCC	1	\$ 1,200.00	\$600.00	AQ2781V High Power OPM & VLS Module	New	Not Discontinued
AQ7280-HE	1	\$ 2,520.00	\$1,550.00	OTDR Main Frame	New	Not Discontinued
AQ2200-331	1	\$ 6,275.00	\$1,255.00	Optical Attenuator Module for AQ2200 Series	Used (Certified Pre-owned)	Discontinued; Mar 31st, 2017;
AQ6370D-22-L1-D/FC/RFC	1	\$ 43,270.00	\$ 33,880.00	AQ6370D Telecom Optical Spectrum Analyzer 600 - 1700 nm (High Performance)	Used (Certified Pre-owned)	Discontinued; Dec 1st, 2023;
Portable Test Instruments						
Model Code	Qty.	List Unit Price	Special Offer Price	Description	Condition	Condition
CA150/TE	4	\$2,244.00	\$1,000.00	HANDY CAL	New	Discontinued; Sept 30th, 2020;
CL420	2	\$840.00	\$170.00	HANDY CAL	Used (Certified Pre-owned)	Not Discontinued